

CYNOSURE®

Palomar  
ConBio

Palomar Icon®

One Platform To Power  
The Future of Aesthetic Treatments



## **Increasing Demand**

Skin revitalization platform revenues are expected to grow by 8.1% each year to an estimated \$478.5 million in 2017<sup>1</sup>



**Before**

Vascular Lesion Removal

**After**

Post 2 Treatments

Khaili Khatri, MD



**Before**

Stretch Marks

**After**

Post 3 Treatments

Sean Doherty, MD



**Before**

Wrinkle Reduction

**After**

Post 1 Treatment

Haneef Alibhai, MD



**Before**

Acne Scars

**After**

Post 3 Treatments

Michael Sinclair, MD



**Before**

Pigmented Lesion Removal

**After**

Post 1 Treatment

Akio Sato, MD

# *Palomar* **ICON**<sup>®</sup>

The next generation of aesthetic light based systems – providing fast treatments with excellent outcomes and user experience.



## **One platform. Multiple applications.**

- Skin Revitalization
- Wrinkles
- Surgical Scars & Acne Scars
- Stretch Marks
- Permanent Hair Reduction
- Vascular & Pigmented Lesions
- Leg Veins

## The Icon system provides the most requested aesthetic treatments with many patented innovations.

- **High Peak Power**  
Optimized results in fewer treatments
- **Advanced Contact Cooling™**  
Increased patient comfort
- **AccuSpectrum™**  
Dual filters minimize unnecessary epidermal exposure for enhanced client comfort
- **Photon Recycling™**  
Reduces energy loss for maximized treatment results
- **Smooth Pulse™**  
Proprietary light delivery system that allows for the usage of greater amounts of light energy
- **Non-Ablative XF Microlens™ & XD Microlens™**  
The XF Microlens provides high non-ablative coverage for fast resurfacing treatments; the XD Microlens offers greater treatment depth for challenging conditions



## The Industries Only Melanin Reader



**Featuring the Skintel® Melanin Reader™** -  
measures melanin content to better  
guide the treatment process.


No more guessing... Skintel provides  
Treatment Confidence



“ “ The 1540 Fractional Laser with XD and XF  
Microlenses has been a wonderful addition to  
my practice. We have seen major improvements  
in acne scars and surgical scars  
in just a few treatments.

-Joel L. Cohen, MD

” ”



An IPL is an essential part of my aesthetic practice making the MaxG my go-to handpiece for unwanted vascular and pigmented lesions. The 1540 fractional handpiece is very effective, reliable and easy to use, making these handpieces a perfect complement on the Palomar Icon aesthetic system.

- Emil Tanghetti, MD



## Palomar Icon Handpieces

	2940	1540	1064+	MaxG™	MaxR™	MaxRs™	MaxYs™
Fractional Skin Resurfacing	●	●					
Permanent Hair Reduction			●		●	●	●
Pigmented & Vascular Lesions and Rosacea				●			
Pigmented Lesions				●			●
Leg Veins			●				
Wrinkles	●						
Melasma		●					
Scars and Stretch Marks		●					




Electrical Requirements	90-264 VAC, 50/60 Hz (20 amp circuit required) Max Power Input 1600 VA
Dimensions	22.9" (58 cm) L x 21.4" (54 cm) H x 40.5" (103 cm) D
Net Weight	135 lbs. (61 kg); modular and mobile design

## One Platform to Power the Future of Aesthetic Treatments

'Skin Rejuvenation: Continually Advancing  
Technology, Medical Insight Inc. 2013: 9.

Please contact your Palomar Icon Consultant  
at 800.886.2966 or visit [www.cynosure.com](http://www.cynosure.com).



### CORPORATE HEADQUARTERS

#### United States

Cynosure, Inc.  
5 Carlisle Road  
Westford, MA 01886 USA  
Tel: 978.256.4200  
Toll-Free: 800.886.2966

[www.cynosure.com](http://www.cynosure.com)

MaxG and MaxY are trademarks and Cynosure, Palomar Icon and Skintel are registered trademarks of Cynosure, Inc.

©2014 Cynosure, Inc. • 921-0485-000 rev 1 03/14